

PRE-LAB ACTIVITY 8

UTILIZING THE WIND TUNNEL SOFTWARE

By: Engr. Arvin N. Gastardo

OBJECTIVES

After completing the calibration, the student will be able to:

1. Utilize the program capabilities.
2. Determine the specific data for different parameters inside the test section.

PROCEDURES

1. Open the program by double clicking the icon  on the desktop, then figure 24 below will be shown.



Figure 24. General Display Model

2. Click the remaining tab pages to find the subject that you want. Each tab will give you an idea about the general description, capabilities, technical data, and the signal diagram of the wind tunnel.
3. Press the lift, drag and pitch button to start the experiment. The experimental screen will show as figure 25 below.




Figure 25. Experiment Screen

Take note:

- No. 1 button: fill in the name of the person using the program
- No. 2 button: key in the room temperature
- No. 3 button: adjust speed control by key in (0-60 Hz) or use by knob button and press. By clicking the start button, the fan will automatically start. Clicking the stop button will stop the fan from rotating.
- No. 4 button: collect the data result
- No. 5 button: Use for deleting the last data
- No. 6 button: Use to delete all recorded data
- No. 7 button: Allows you to select any preset test subject.
- No.8 button: Will allow you to save all logged data. The program will show “Data sheet” that its format is as figure 3 below.
- No. 9 button: Click to present analysis data on screen and include this data in the saved data file

Experiment : Lift , Drag and Pitching						
Tested by : Date : 06/09/56 Time : 10:41						
Model : Aerofoil , Dimensions: 100 x 100 mm						
Air Velocity 24.5 m/s Room Temperature: 30.0 Degree Celsius						
Angle of Attack (degree)	Lift (N)	Drag (N)	Pitching (N.m)	Coefficient		
				CL	CD	Cm
1	0.10	0.16	-0.031			
10	1.34	0.47	-0.025			
20	2.06	1.15	0.002			
30	2.33	1.72	-0.018			
40	2.52	2.34	-0.031			
50	2.46	3.08	-0.034			
59	2.06	3.70	-0.041			
70	2.10	4.38	-0.057			
80	1.25	4.79	-0.060			
90	0.07	4.76	-0.054			

Figure 26. Data Results

- 4. Key in your full name.
- 5. Key in the room temperature.

6. Select model for the experiment according to the experiment setup. We will use the long aerofoil for this experiment.
7. Adjust the wind velocity (20 m/s) control for the experiment.
8. Set the model's angle of attack to zero degree
9. Collect the data result
10. Rotate the model at 10 degrees increment and repeat step 9.
11. Click the save to file button and choose the desktop as file location. Assign your name as the file name. Provide data sheet on the next page under table A-1.
12. Write your conclusion for this activity on the last page.

PRE-LAB ACTIVITY 2 RESULTS SHEET

NAME OF STUDENT	
DATE	
COURSE NUMBER/COURSE NAME	
SEMESTER & ACADEMIC YEAR	
NAME OF INSTRUCTOR/PROFESSOR	

Table A-1

Patient Information	
First Name	
Last Name	
Room Number	
Ward	
Age	
Gender	
Referral Source	
History of Present Illness	
Onset of symptoms	
Duration of symptoms	
Frequency of symptoms	
Severity of symptoms	
Associated symptoms	
Previous medical history	
Family history	
Social history	
Review of Systems	
General	
Cardiovascular	
Respiratory	
Gastrointestinal	
Genitourinary	
Neurological	
Musculoskeletal	
Dermatological	
Endocrine	
Hematological	
Immunological	
Psychological	
Other	
Physical Examination	
Vital Signs	
General Appearance	
Head and Neck	
Chest	
Abdomen	
Genitourinary	
Neurological	
Musculoskeletal	
Dermatological	
Endocrine	
Hematological	
Immunological	
Psychological	
Other	
Laboratory and Diagnostic Tests	
Blood Tests	
Urine Tests	
Imaging Studies	
Other Tests	
Treatment and Management	
Medications	
Surgery	
Physical Therapy	
Other Treatments	
Follow-up and Prognosis	
Next Appointment	
Patient Education	
Discharge Instructions	
Referral to Specialist	
Other Notes	

CONCLUSION: